

Software Quick Start (SDK)



1



Insert the **EyeTech USB Drive** into an available **USB Port** on your computer.

2

Double-click on the Installer Icon on the flash drive and follow the on-screen prompts. Accept all defaults. The installer will place the SDK files in: C:\Program Files\EyeTechDS\QuickLINK_x.x.x.x by default.

3

To verify camera operation, run the Quick Glance.exe program in:
C:\Program Files\EyeTechDS\QuickLINK_x.x.x.x\bin

For additional help, visit our Question and Answer website: <http://answers.eyetechds.com/>